



The DualFlex - when Comfort is a priority!!

DualFlex is a HYBRID mattress replacement system, combining both high performance static foam, and active alternating air cell technology for exceptional patient comfort and healing.

OCTOBER SALE - 10% off Dual Flex - pumps and mattresses can be purchased separately.

[Contact us](#) for pricing and further product queries.



KEY FEATURES:

- * **Made in the UK**
- * **SWL** - up to 250kgs
- * **Risk** - High
- * **Ultra Quiet Pump** - less than 28db, for a peaceful sleeping environment.
- * **Auto-set** - no need to input weight, with ability to manually adjust for patient comfort
- * **Patient Transfer mode**
- * **GELTEX** inside represents the next generation of medical mattress technology. Its patented formula provides a unique combination of breathability, pressure distribution and body support. **GELTEX** is created by fusing elastic gel with air-permeable foam, resulting in a breathable, open cell material that promotes and improved microclimate between patient and mattress."
<https://recticelflexiblefoams.com/geltex>
- * **Nylon Air cells** with TPU lamination for longevity
- * **Low Interface Pressures** for improved blood circulation and wound healing.
- * **Multi Zoned** - Head, sacrum and heel zones
- * **Channels for hose and cell connections** to avoid pressure points on the surfaces.
- * **Flexible contour design** for profiling beds
- * **IPS Technology** - perfectly adjusts air pressure to suit individual patients and their movements.
- * **Handles** for improved infection control
- * **Hidden zip** for safely securing it to the bed frame.
- * **Buckle and straps** for ease of carrying.

