

# 7

# carilex<sup>®</sup>

## Low Profile - Alternating Mattress Replacement



Mattress Replacement



Durable Material



Multiple Mattress Sizes



Cable Management



Easy Cleaning



5 inch Air Cell



Protective PE Layer

Global Leader in Wound Care, Distributed in Over 40 Countries Worldwide

# 7 Low Profile - Alternating Mattress Replacement

The Carilex Mattress Replacements are designed to facilitate blood circulation and provide maximum tissue pressure redistribution and relief for treating and preventing pressure injuries.

They form the ideal therapy system when used with a compatible pump providing low noise levels and the comfortable modular cell design for superior sleeping quality.

## GENERAL FEATURES

Patient Envelopment	Low profile, alternating air mattress replacement - envelops and off loads patient body mass by expanding the therapy contact area and depth
Quick Connect Coupling	Easy connection between power unit and mattress
Cable Management	Positioned on one side for no fuss cable and cord management
CPR	Pull Tab deflates the entire mattress within seconds
Washable Cover	2 way stretch
Zip on Foam Base	Extra 2" foam protection for patient support

## OPTIONAL EXTRA

Four Direction Stretch Coverlet	High quality DARTEX washable four-way stretch, shear force reduction, excellent moisture permeability
---------------------------------	---

## TECHNICAL DETAILS

Risk/Category	High
Waterlow Score	15+
Pressure Sore Grading	2
Safe Working Load	Up to 150kg
Air Cell Configuration	16 x 5" cells
Cycle Times	10 minutes, 2:1
Available Size/Dimensions	Single 90cm x 200cm x 18cm, Euro KS 100cm x 200cm x 18cm, King Single 107cm x 200cm x 18cm
Securing Mattress to Bed	Buckle straps
Travel/Weight	Mattress rolls up into a bag, weighing 4.3kg
Clinical Settings	Suitable for Hospital, Home Care and Nursing Home environments
Compatible Pumps	Entrix NX and Centrius
Warranty	2 year

# carilex<sup>®</sup>

Manufacturing for over 30 Years - 10 Years in the Australian Market