

Centrius

carilex®

Pump

The Intelligent Pressure Sensing Technology (IPS), responds to patient movements on the mattress by automatically adjusting internal mattress pressure, this allows Centrius to regulate airflow and continuously provide total envelopment.



Alternating
Pressure Relief



MAX Inflation



Whisper Quiet



CLP Therapy
Pressure Redistribution



Lightweight Power Unit



Memory Recall

Global Leader in Wound Care, Distributed in Over 40 Countries Worldwide



GENERAL FEATURES

1 Power/Mute/Lock	
2 Panel Lock	System is locked when this light is orange
3 Comfort Level Settings	Auto adjust to individual patient comfort
4 Mode Button	Changes the therapy mode
5 CLP Therapy	All air cells maintain constant pressure support
6 Max Inflate	Rapidly hyper-inflates the mattress providing a firm and convenient surface for nursing tasks
7 Dynamic Mode	2:1 alternating therapy, 10 minute cycle
8 Low Pressure Alert	Orange LED indicates insufficient pressure
Intelligent Pressure Sensing	Responds to patient movements on mattress by automatically adjusting internal mattress pressure
Power Failure Alarm	Audible and visual power failure alarms to alert caregivers in case of abnormal power outage
Hands Free Automatic System	Soft, Medium or Firm fully automated comfort modes, with auto upright support
Transport Mode	Connect the two quick couplings together to retain air in the mattress, can stay inflated for up to 15 hours

TECHNICAL DETAILS

Risk/Category	High to Very High
Waterlow Score	15+ to 20+
Pressure Sore Grading	2 to 4
Safe Working Load	0 to 250kg
Compatible Support Surfaces	5" Overlay, 7", Lux and FIT Mattress Replacements - Single and King Single
Travel/Weight	Pump weighs 2.3kg
Clinical Settings	Suitable for Hospital, Home Care and Nursing Home environments
Warranty	2 year

carilex[®]

Manufacturing for over 30 Years - 10 Years in the Australian Market