

# Coziny 100 <sup>®</sup> carilex

Paediatric Pressure Care Support Surface



Overlay



1 inch Air Cells



Modular Air Cell  
Construction



Quick Connect  
Coupling



DTX Coverlet



Alternating  
Pressure Relief

Global Leader in Wound Care, Distributed in Over 40 Countries Worldwide

Coziny 100 mattresses are made entirely with Extra Soft 100% TPU material to ensure baby's safety and comfort. The Extra Soft 100% TPU fabric air inlets also protect baby from sharp edges and protruding objects.

### GENERAL FEATURES

Extra Soft 100% TPU Material	Ensuring patients safety from sharp edges or protruding objects while providing optimal comfort
Special Temperature Distributing Air Channels	Allowing the heat from the incubator to evenly keep the baby warm throughout
Hidden Zipper	With protection pad for baby's safety and better comfort
Overlay	To be placed on top of incubator mattress
Weight Capacity	For patients up to 3.4kgs
Dartex Washable Cover	4 way stretch

### TECHNICAL DETAILS

Risk/Category	Medium to High
Waterlow Score	15+
Pressure Sore Grading	3
Safe Working Load	Up to 3.4kg
Air Cell Configuration	40 x 1" cells
Cycle Times	5, 10 or 15 minutes
Available Size/Dimensions	Incubator 65cm x 35cm x 3cm
Travel/Weight	0.5kg
Clinical Settings	Suitable for Hospital and Home Care environments
Compatible Pumps	Coziny Plus
Warranty	2 year

carilex<sup>®</sup>

Manufacturing for over 30 Years - 10 Years in the Australian Market