

Coziny 300 [®] carilex

Paediatric Pressure Care Support Surface



Mattress Replacement



Modular Air Cell Construction



DTX Coverlet



1.5 inch Air Cells



Quick Connect Coupling



Alternating Pressure Relief

Global Leader in Wound Care, Distributed in Over 40 Countries Worldwide

Coziny Plus offers multiple types of mattresses designed with extra attention to details, providing remarkable comfort. Extra consideration was taken to eliminate sharp or protruding components from the mattresses, ensuring that we deliver the most protective wound care prevention and therapy for the most sensitive skin.

GENERAL FEATURES

Three Connecting Hose Outlets.	The three positions where the connecting hoses can be attached provides compatibility with a wide range of baby cribs
1.5cm Foam Base	Integrated for extra protection in case of power outage
Weight Capacity	For patients up to 25kgs
Dartex Washable Cover	4 way stretch

TECHNICAL DETAILS

Risk/Category	Medium to High
Waterlow Score	15+
Pressure Sore Grading	3
Safe Working Load	Up to 25kg
Air Cell Configuration	29 x 1.5" cells
Cycle Times	5, 10 or 15 minutes
Available Size/Dimensions	Baby Cot 140cm x 70cm x 8cm
Travel/Weight	3.5kg
Clinical Settings	Suitable for Hospital and Home Care environments
Compatible Pumps	Coziny Plus
Warranty	2 year

carilex[®]

Manufacturing for over 30 Years - 10 Years in the Australian Market