

FIT

carilex[®]

Deep Cell - Alternating Mattress Replacement



Mattress Replacement



Durable Material



Multiple Mattress Sizes



Cable Management



Easy Cleaning



8 inch Air Cell



Protective PE Layer

Global Leader in Wound Care, Distributed in Over 40 Countries Worldwide

FIT Deep Cell - Alternating Mattress Replacement

The Carilex Mattress Replacements are designed to facilitate blood circulation and provide maximum tissue pressure redistribution and relief for treating and preventing pressure injuries.

They form the ideal therapy system when used with a compatible pump providing low noise levels and the comfortable modular cell design for superior sleeping quality.

GENERAL FEATURES

Patient Envelopment	Low profile, alternating air mattress replacement - envelops and off loads patient body mass by expanding the therapy contact area and depth
Quick Connect Coupling	Easy connection between power unit and mattress
Cable Management	Positioned both sides, for no fuss cable and cord management
CPR	CPR Dial deflates the entire mattress within seconds
Washable Cover	4 way stretch - high quality DARTEX washable four-way stretch, shear force reduction, excellent moisture permeability

OPTIONAL EXTRAS

Zip on Foam Base	Extra 2" foam protection for patient support
------------------	--

TECHNICAL DETAILS

Risk/Category	Very High
Waterlow Score	20+
Pressure Sore Grading	4
Safe Working Load	Up to 250kg
Air Cell Configuration	20 x 8" square cells
Cycle Times	10 minutes, 2:1
Available Size/Dimensions	Single 90cm x 200cm x 20cm King Single 107cm x 200cm x 20cm
Securing Mattress to Bed	Buckle straps
Travel/Weight	Mattress rolls up into a bag, weighing 9kg
Clinical Settings	Suitable for Hospital, Home Care and Nursing Home environments
Compatible Pumps	Centrius
Warranty	2 year

carilex[®]

Manufacturing for over 30 Years - 10 Years in the Australian Market