

LUX

carilex®

Low Profile - Alternating Mattress Replacement



Mattress Replacement



Durable Material



Multiple Mattress Sizes



Cable Management



Easy Cleaning



7 inch Air Cell



Protective PE Layer

Global Leader in Wound Care, Distributed in Over 40 Countries Worldwide

LUX Low Profile - Alternating Mattress Replacement

The Carilex Mattress Replacements are designed to facilitate blood circulation and provide maximum tissue pressure redistribution and relief for treating and preventing pressure injuries.

They form the ideal therapy system when used with a compatible pump providing low noise levels and the comfortable modular cell design for superior sleeping quality.

GENERAL FEATURES

| | |
|------------------------|--|
| Patient Envelopment | Low profile, alternating air mattress replacement - envelops and off loads patient body mass by expanding the therapy contact area and depth |
| Quick Connect Coupling | Easy connection between power unit and mattress |
| Cable Management | Positioned side for no fuss cable and cord management |
| CPR | CPR Dial deflates the entire mattress within seconds |
| Washable Cover | 2 way stretch |

OPTIONAL EXTRAS

| | |
|---------------------------------|---|
| Zip on Foam Base | Extra 2" foam protection for patient support |
| Four Direction Stretch Coverlet | High quality DARTEX washable four-way stretch, shear force reduction, excellent moisture permeability |

TECHNICAL DETAILS

| | |
|---------------------------|--|
| Risk/Category | Very High |
| Waterlow Score | 20+ |
| Pressure Sore Grading | 4 |
| Safe Working Load | Up to 180kg |
| Air Cell Configuration | 20 x 7" square cells |
| Cycle Times | 10 minutes, 2:1 |
| Available Size/Dimensions | Single 90cm x 200cm x 18cm King Single 107cm x 200cm x 18cm |
| Securing Mattress to Bed | Buckle straps |
| Travel/Weight | Mattress rolls up into a bag, weighing 7kg |
| Clinical Settings | Suitable for Hospital, Home Care and Nursing Home environments |
| Compatible Pumps | Entrix NX and Centrius |
| Warranty | 2 year |

carilex[®]

Manufacturing for over 30 Years - 10 Years in the Australian Market