

MobiCare Plus

Wheelchair or Chair Cushion



GENERAL FEATURES

Extra Support at the Sacrum Area	Specially designed deep air cells for excellent immersion and pressure relief.
Durable Material	High quality material intended for long term usage.
Quick Release Coupling	Provides an easy connection between power unit and cushion
Better Immersion	10 x 3.5" deep air cells

TECHNICAL DETAILS

Risk/Category	Med - High
Waterlow Score	15+
Pressure Sore Grading	3
Safe Working Load	Up to 150kgs
Cushion Dimensions	36cm x 45cm x 9cm, 45cm x 45cm x 9cm or 60cm x 45cm x 9cm
Travel/Weight	Weighing 1.5kg
Clinical Settings	Suitable for Hospital and Home Care environments
Compatible Pumps	CariChair
Warranty	1 year